


Extended analysis of a validated urine-exosome signature to predict high grade prostate cancer on initial biopsy maintains performance across multiple sub-groups.

Michael Donovan (1), Stefan Bentink (2), Mikkel Noerholm (2), Vince O'Neill (3), Johan Skog (3);

(1) Icahn School of Medicine at Mount Sinai, New York, NY; (2) Exosome Diagnostics GmbH, Martinsried, Germany, (3) Exosome Diagnostics, Cambridge, MA

Introduction: We recently completed a prospective observational clinical trial of a first-catch urine exosome gene expression signature ExoIntelliScore Prostate (formerly known as EXO106) to predict high-grade prostate cancer, HGPCA (> / = GS7) on initial biopsy. The assay accurately predicted HGPCA with an NPV > 90% for an intended use population: men > / = 50 years with an equivocal PSA of 2-10 ng/mL presenting for their first biopsy. Given the current high prostate biopsy rate and potential impact of race we sought to further understand performance of the assay by constructing sub-cohorts from the validation study.



Methods: Utilizing the ExoIntelliScore Prostate result and clinical data from the 519 patient intended use clinical validation cohort (IU) we constructed several sub-groups including: i. prior negative (A, n = 149), ii. combined prior negative + initial biopsy (B, n = 668) and African Americans, AA (C, n = 87). Area under the curve (AUC), NPV, PPV, sensitivity and specificity with validated cut-point assess performance.

Results: There was demographic comparability between the initial intended use (IU; Figure 1) and various subgroups, i.e. prior negative (A; Figure 2), combined prior negative and initial biopsy (B; Figure 3) and African Americans (C; Figure 4) with some exceptions observed (bold) for AA patients; family history (IU, 23%; A: 27%, B: 24%, C: **30%**; positive biopsy rate (IU, 48%; A: 34%, B: 45%, **C: 52%**; and > / = GS7 (IU, 28%; A: 13%, B: 27%; C: **34%**). The AUC range for ExoIntelliscore Prostate + standard of care (i.e. age, race, PSA and family history) between IU and A, B subgroups was 0.72-0.74; AUC for AA was 0.62, except for prior negative biopsy AA patients, AUC 0.80; supporting importance of prior negative biopsy. Of note, applying the ExoIntelliscore Prostate cut-point, A, B yielded an NPV of 91%; including NPV of 89% for AA (initial biopsy) and NPV of 91% for AA with initial + prior negative biopsy.

Figure 1. Intended Use (IU) Initial biopsy cohort, serum PSA 2-10 ng/mL; n=519
Discrimination for High Grade Cancer; SOC=Age, PSA, Fam History, Race)




Figure 3. Group B: Initial + First Biopsy, serum PSA 2-10ng/mL; n=668




Figure 2. Group A: Repeat Biopsy, serum PSA 2-10ng/mL; n=149




Figure 4. Group C: African American, initial biopsy, serum PSA 2-10ng/mL; n=87




Table 1. Comparison of Intended use cohort with various sub-groups.

Cohorts	N	HGP	TP	TN	FP	FN	NPV	PPV	SENS	SPEC
IU - Initial biopsy, (PSA 2-10)	519	28.5	136	126	245	12	91.3	35.7	91.9	34.0
A -Repeat Biopsy; (PSA 2-10)	149	12.7	16	33	97	3	91.7	14.2	84.2	25.4
B - First and Repeat Biopsy (PSA 2-10)	668	25.0	152	155	346	15	91.2	30.5	91.0	31.0
C - African-American; Initial Biopsy (PSA 2-10)	87	33.3	27	16	42	2	88.9	39.1	93.1	27.6

HGP, high Gleason grade (>/-GS7) prevalence; TP, true positive; TN, true negative; FP, false positive, FN, false negative; NPV, negative predictive value; PPV, positive predictive value; SENS, sensitivity; SPEC, specificity.

Conclusions: The ExoIntelliScore Prostate performed equally well in men with or without a prior negative biopsy with comparable results also seen for African Americans. Additional confirmatory studies with more patients are necessary to confirm initial observations

References:

- Miranda KC, Bond DT, McKee M, Skog J, Paunescu TG, Da Silva N et al. Nucleic acids within urinary exosomes/microvesicles are potential biomarkers for renal disease. *Kidney Int.* 2010;78:191-199.
- Donovan MJ, Noerholm M, Bentink S, Belzer S, Skog J et al. A molecular signature of PCA3 and ERG exosomal RNA from non-DRE urine is predictive of initial prostate biopsy result. *Prostatic Diseases and Prostate Cancer.* 2015, in press
- He J, Schepmoes AA, Shi T, Wu C, Fillmore TL, Gao Y et al. Analytical platform evaluation for quantification of ERG in prostate cancer using proteins and mRNA detection methods. *J Transl Med.* 2015;54:1-14.
- Klotz L. Active Surveillance for Low-Risk Prostate Cancer. *Curr Urol Rep* 2015; 16: 1-10.
- Nilsson J, Skog J, Norstrand A, Baranov B, Mincheva-Nilsson L, Breakefield XO, Widmark A; Prostate Cancer-Derived Urine Exosomes; A Novel Approach to biomarkers for Prostate Cancer; *British Journal of Cancer* (2009) 100, 1603-1607

